

schülke -+



Ready-to-use alcohol-based rapid disinfection for non-invasive medical products and other surfaces.

pursept® A Xpress

Our Plus

- ready to use and highly effective in only 15 seconds
- broad biocidal activity
- dries rapidly with hardly leaving any smears
- very fresh fragrance
- dermatologically tested

Application areas

Alcohol-based rapid disinfection of medical devices and other surfaces in all areas with an increased risk of infection and where short contact times are required, e.g.:

- Patient nearby surfaces
- Patient treatment units
- Examination couches
- Operating tables and nearby work surfaces
- Surfaces of medical equipment



Safety in
a few seconds.

Microbiological efficacy

Efficacy	Concentration	Contact time
bactericidal <i>in accordance with VAH Guidelines with short contact times - clean conditions</i>	ready-to-use	15 sec.
tuberculocidal <i>in accordance with VAH Guidelines with short contact times - dirty conditions</i>	ready-to-use	30 sec.

Efficacy	Concentration	Contact time
levurocidal <i>in accordance with VAH Guidelines with short contact times - clean conditions</i>	ready-to-use	15 sec.
enveloped viruses (incl. HIV, HBV, HCV) <i>in accordance with DVV (German Registered Association for Combating Viral Diseases)//RKI Guideline</i>	ready-to-use	15 sec.
Adenovirus <i>in accordance with DVV (German Registered Association for Combating Viral Diseases)//RKI Guideline</i>	ready-to-use	1 min.
Norovirus <i>EN14476 - dirty conditions</i>	ready-to-use	15 sec.
Rotavirus <i>in accordance with DVV (German Registered Association for Combating Viral Diseases)//RKI Guideline</i>	ready-to-use	15 sec.

Certificates

- IHO listed
- VAH certificate



CE 0297



Product data

Composition:

100 g agent solution contain the following active ingredients: 58,6 g Ethanol (94%) , 0,03 g Amines, N-C12-C14 (even-numbered)-alkyltrimethylenedi-, reaction products with chloroacetic acid.

Labelling according to Regulation (EC) No 648/2004: perfume

Chemical-physical data

Color	colourless
Density	ca. 0,90 g/cm ³ / 20 °C
Flash point	23 °C
Form	liquid
pH	7,5 / 20 °C / (undiluted)
Viscosity, dynamic	No data available

Special advice

Instruction for use

pursept® A Xpress is a ready-to-use solution. Apply the undiluted pursept® A Xpress onto objects and surfaces, wipe and allow it to take effect. Always aim for complete surface coverage. Maximum application rate 50 ml/m². Make sure that all visible soiling is removed before disinfection. Do not treat surfaces that are sensitive to alcohol, e.g. acrylic glass. Not for final disinfection of semicritical and critical medical devices!

Use disinfectants safely. Always read the label and product information before use.

The product is highly compatible with metals and synthetics (except acrylic glass and alcohol-sensitive varnishes).

Wearing gloves is recommended. Do not use near ignition sources.

Please follow carefully any national fire and explosion protection guidelines when using alcohol-based disinfectants.

Information for order

Item	Delivery form	Item no.
pursept® A Xpress 1 l spray bottle	10/ Carton	on request
pursept® A Xpress 5 l	1/ Canister	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

Accessories

Accessories	Item no.
Bed bracket for 1 l square bottle (curved)	134415
Bed bracket for square bottle 1 l (straight)	134413
Can key for 5 + 10 l	135810
Decanter funnel	117901
measuring cup 500 ml	136101
Opener all-in-one (1l eurobottle/can/container)	70000804
schülke dosing feeder 5 l / 10 l (20ml / stroke)	117101
schülke tap for 5 l / 10 l can	135501
wall bracket plastics 1l	134416

Related Products

- pursept® A Xpress wipes

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining out high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: www.schuelke.com.

For individual questions:

Customer Sales Service

Phone: +49 40 52100-666

E-Mail: info@schuelke.com



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.



schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 (0) 40 - 52100 - 0
Fax +49 (0) 40 - 52100 - 318
www.schuelke.com
mail@schuelke.com