LiquiSampler | sample collector for liquids











LiquiSampler for liquids, according to DIN 53 242. Made of ultra pure and chemically inert PTFE/FEP or PP for contamination-free sampling.
Sampling from open and closed containers, barrels, tanks, silos, water courses. Suitable for all commercially available barrels and containers with opening of at least 25 mm in diameter.

Push-buttons colour-coded to prevent mix-up, available as optional accessories.
Easy cleaning - all surfaces are free of pores and crevices, preventing accumulation of dirt. Only round screw threads proven in foodstuffs hygiene are used. Disassembly and cleaning are extremely simple.

- Practical one-handed use at the push of a button
- With scraping ring for marking of depth and cleaning
- Transparent
- Available in PTFE/ FEP or PP transparent
- Point sampling, all-layer and soil samples possible

Use

All-layer sample

1. Slowly insert the LiquiSampler into the liquid with the button pressed.
2. Release the push-button, the valve closes.
3. Withdraw the lance and take out the sample.

Point sampling/soil sample

1. Insert the LiquiSampler to the targeted point/to the ground with closed valve.
2. Press the button - the liquid flows in.
3. Release the push-button, the valve closes.
4. Pull out the LiquiSampler, take out the sample.

LiquiSampler

Material	Length cm	Dia. mm	Volume ml	Price €	Item no.
PP, transparent	60	25	150	139.35	5330-1060 N4
PP, transparent	100	25	250	185.45	5330-1100 N4
PP, transparent	200	25	500	277.50	5330-1200 N4
PTFE/FEP, transparent	60	25	150	215.35	5330-2060 N4
PTFE/FEP, transparent	100	25	250	287.75	5330-2100 N4
PTFE/FEP, transparent	200	25	500	397.10	5330-2200 N4

Accessories

Description	Price €	Item no.
Push-buttons for set LiquiSampler in blue, green, yellow, white and black	24.05	5330-0149 N4
Cleaning brushes 25 mm dia., length 100 cm	35.60	5304-0100 N4