ViscoSampler | sample collector for medium viscosity substances









ViscoSampler with powerful suction for viscous media. The media is drawn up and discharged using the suction flask. Made of ultra-pure, chemically inert PTFE/FEP for contamination-free sampling.

Specially developed for viscous substances such as slurries, silt, wet clay and soil samples, soil sediments, oils, emulsions, creams etc.

Colour-coded handles to prevent mix-up, available as optional accessories.

Easy cleaning - all surfaces are poreless and without crevices. Dirt accumulation is prevented. Only round screw threads tried and tested in the field of food hygiene are used. Disassembly and cleaning are extremely simple.



QualiSampler set
The handy set, which comes complete with carrying case, combines the LiquiSampler and the ViscoSampler with practical accessories.
Models in PP or PTFE/ FEP.

- QualiSampler set consisting of:
- LiquiSampler, 100 cm
- ViscoSampler, 100 cm
- Case, five bottles, cleaning brush

ViscoSampler

Material	Length cm	Dia. mm	Capacity ml	Price €	Item no.
PP, transparent	60	25	160	215.35	5331-1060 N4
PP, transparent	100	25	300	254.35	5331-1100 N4
PP, transparent	200	25	650	377.55	5331-1200 N4
PTFE/FEP, transparent	60	25	160	288.85	5331-2060 N4
PTFE/FEP, transparent	100	25	300	328.00	5331-2100 N4
PTFE/FEP, transparent	200	25	650	463.75	5331-2200 N4

Accessories

Description	Price €	Item no.
Handles ViscoSampler set in blue, green, yellow, white and black	39.85	5331-0020 N4
Cleaning brushes 25 mm dia., length 100 cm	35.60	5304-0100 N4
Drain safety stand, see page 55	42.25	5618-1000 N4

QualiSampler set

Material	Price €	Item no.
PP	471.95	5330-1000 N4
PTFE / FEP	662.95	5330-1001 N4

LiquiSampler | sample collector for liquids











LiquiSampler for liquids, according to DIN 53 242. Made of ultra pure and chemically inert PTFE/FEP or PP for contamination-free sampling.
Sampling from open and closed containers, barrels, tanks, silos, water courses. Suitable for all commercially available barrels and containers with opening of at least 25 mm in diameter.

Push-buttons colour-coded to prevent mix-up, available as optional accessories.
Easy cleaning - all surfaces are free of pores and crevices, preventing accumulation of dirt. Only round screw threads proven in foodstuffs hygiene are used. Disassembly and cleaning are extremely simple.

- Practical one-handed use at the push of a button
- With scraping ring for marking of depth and cleaning
- Transparent
- Available in PTFE/ FEP or PP transparent
- Point sampling, all-layer and soil samples possible

Use

All-layer sample

1. Slowly insert the LiquiSampler into the liquid with the button pressed.
2. Release the push-button, the valve closes.
3. Withdraw the lance and take out the sample.

Point sampling/soil sample

1. Insert the LiquiSampler to the targeted point/to the ground with closed valve.
2. Press the button - the liquid flows in.
3. Release the push-button, the valve closes.
4. Pull out the LiquiSampler, take out the sample.

LiquiSampler

Material	Length cm	Dia. mm	Volume ml	Price €	Item no.	
PP, transparent	60	25	150	139.35	5330-1060	N4
PP, transparent	100	25	250	185.45	5330-1100	N4
PP, transparent	200	25	500	277.50	5330-1200	N4
PTFE/FEP, transparent	60	25	150	215.35	5330-2060	N4
PTFE/FEP, transparent	100	25	250	287.75	5330-2100	N4
PTFE/FEP, transparent	200	25	500	397.10	5330-2200	N4

Accessories

Description	Price €	Item no.	
Push-buttons for set LiquiSampler in blue, green, yellow, white and black	24.05	5330-0149	N4
Cleaning brushes 25 mm dia., length 100 cm	35.60	5304-0100	N4