

Metal caps

For culture flasks and culture medium bottles.

DURAN Group

For neck diam. mm	Material	PK	Cat. No.
38	Stainless steel	1	9.010 481 1
38	Aluminium, anodised blue	1	9.010 482 2



9.010 481

9.010 482

3 Metal caps, Erlenmeyer, Aluminium

Anodised in a choice of colours.

schuett-biotec

Special caps for sterile, but not hermetically sealing, closure of culture flasks and media bottles with straight, rimmed or rimless necks as indicated. With nichrome steel springs that clamp onto the outer wall of the container.

Cleaning advice available on request.

For neck diam. mm	Colour	With rim	PK	Cat. No.
37/39	silver	no	10	9.231 456
37/39	blue	no	10	9.231 457
37/39	red	no	10	9.231 458



3

4 Metal caps, Labocap without handle

schuett-biotec

For neck diam. mm	Colour	Handle	PK	Cat. No.
9/10	blue	without	100	9.231 421
12/13	blue	without	100	9.231 422
14/15	blue	without	100	9.231 423
15/16	blue	without	100	9.231 424
17/18	blue	without	100	9.231 425
19/20	blue	without	100	9.231 426
21/23	blue	without	100	9.231 427
24/26	blue	without	100	9.231 428
28/30	blue	without	100	9.231 429
9/10	red	without	100	9.231 441
12/13	red	without	100	9.231 442
14/15	red	without	100	9.231 443
15/16	red	without	100	9.231 444
17/18	red	without	100	9.231 445
19/20	red	without	100	9.231 446
21/23	red	without	100	9.231 447
24/26	red	without	100	9.231 448
28/30	red	without	100	9.231 449
12/13	yellow	without	100	9.231 460
15/16	yellow	without	100	9.231 461
17/18	yellow	without	100	9.231 462
19/20	yellow	without	100	9.231 339
12/13	green	without	100	9.231 463
15/16	green	without	100	9.231 464
17/18	green	without	100	9.231 465
19/20	green	without	100	9.231 340
09/10	silver	without	100	9.231 369
12/13	silver	without	100	9.231 372
14/15	silver	without	100	9.231 374
15/16	silver	without	100	9.231 375
17/18	silver	without	100	9.231 377
19/20	silver	without	100	9.231 379
21/23	silver	without	100	9.231 381
24/26	silver	without	100	9.231 385
28/30	silver	without	100	9.231 388
12/13	black	without	100	9.231 334
15/16	black	without	100	9.231 466
17/18	black	without	100	9.231 337
19/20	black	without	100	9.231 341
12/13	violet	without	100	9.231 335
15/16	violet	without	100	9.231 336
17/18	violet	without	100	9.231 338
19/20	violet	without	100	9.231 342



4



1 Metal caps, Labocap with handle, Aluminium

Anodised in a choice of colours. For sterile, but not hermetically sealing, closure of containers. With nichrome steel springs that clamp onto the outer wall of the container. Autoclavable. Cleaning advice available on request.

schuett-biotec

For neck diam. mm	Colour	Handle	PK	Cat. No.
15/16	silver	with	100	9.231 345
17/18	silver	with	100	9.231 347
15/16	blue	with	100	9.231 403
17/18	blue	with	100	9.231 404
15/16	red	with	100	9.231 413
17/18	red	with	100	9.231 414
15/16	yellow	with	100	9.231 326
17/18	yellow	with	100	9.231 330
15/16	green	with	100	9.231 327
17/18	green	with	100	9.231 331
15/16	black	with	100	9.231 328
17/18	black	with	100	9.231 332
15/16	violet	with	100	9.231 329
17/18	violet	with	100	9.231 333



2 Aluminium caps, pure aluminium

Matt finish. Without internal clips. Loose-fit caps for sterile, but not hermetically sealing, closure of culture flasks and other containers. Can also be used as open specimen holders. Autoclavable.

schuett-biotec

Int. diam. mm	Height mm	PK	Cat. No.
11	25	100	9.230 821
12	25	100	9.230 822
13	25	100	9.230 823
14	25	100	9.230 824
15	25	100	9.230 825
16	30	100	9.230 826
17	30	100	9.230 827
18	30	100	9.230 828
20	30	100	9.230 829
22	30	100	9.230 830
28	30	100	9.230 831
32	30	100	9.230 832
39	30	10	9.230 833
44	40	10	9.230 834
48	40	10	9.230 835
55	40	10	9.230 836
58	40	10	9.230 837

We can **supply** this
manufacturer's
whole
product range !



*schuett***biotec**.de