

Shakers, mixers and stirrers » MMS-3000 magnetic stirrer

MMS-3000 mini magnetic stirrer

Compact magnetic stirrer for routine laboratory procedures.

- Adjustable speed control: 0 to 3000 rpm
- Maximum stirring volume: 20 litres

Handles volumes up to 20 litres

Rugged stainless steel plate provides a working surface of 160 mm

Quiet operation, small size – fits neatly into the workspace

MMS-3000 is equipped with a detachable stand that allows inserting different sensors (temperature, pH etc.) inside the liquid.

Magnetic cylinder-shape stirring bar (6 x 25 mm) supplied as standard

Easy to use graduated dial for quick and convenient selection of mixing speed



- Applications:**
- PH metering, extraction and dialysing with small quantities of substances
 - Any requirement for stirring solutions up to 20 litres

Shakers, mixers and stirrers » MMS-3000 » Models and specifications

Stirrers – models and specifications

		MMS-3000	
		Mini magnetic stirrer	
		 h: 75 mm d: 230 mm w: 185 mm weight: 1.5 kg	
Ambient temperature range	°C	+4 to 40	
Speed	rpm	0 to 3000	
Liquid stirring viscosity	mPa	up to 1170	
Stirring volume	L	20	
Plate dimensions	mm	ø160	
Input voltage	V dc	12	
Input current	A	0.3	