

Vertical aluminium autoclaves

CertoClav CV EL



Autoclave

8501312(12L) 8501312A(12L)
8501318(18L) 8501318A(18L)

Description:

- Sterilisation up to 140°C, proven design since 50 years
- Flow evacuation procedure
- Process control by thermostat and precision valve
- Precision thermometer
- Precision pressure gauge
- Aluminium chamber
- Available in 12 and 18 litre volume
- Max. temperature: 140°C
- Max. pressure (rel.): 2.7 bar (39 psi)
- 220-240Vac, 50Hz, 1900W
- Certificate: CE, TÜV

Included:

Autoclave, base plate, mains cable

Weight and dimensions; model specific differences:

Model	Sterilisation temperature	Dimension chamber DxD	Dimension package WxHxD	Weight package
CV EL 12L 125°C/140°C	Adjustable to 125°C or 140°C	24 x 26 cm	63 x 40 x 37 cm	8 kg
CV EL 12L 115°C/121°C	Adjustable to 115°C or 121°C	24 x 26 cm	63 x 40 x 37 cm	8kg
CV EL 18L 125°C/140°C	Adjustable to 125°C or 140°C	24 x 39 cm	63 x 40 x 37 cm	9kg
CV EL 18L 115°C/121°C	Adjustable to 115°C or 121°C	24 x 39 cm	63 x 40 x 37 cm	9kg



Wire basket
8583010



Condensation vessel
(5L) with hose (1,0m)
8582020

Type	Article name	Art. No.		
Autoclave	CertoClav EL 12L GS 125/140°C	8501312		
Autoclave	CertoClav EL 12L GS 115/121°C	8501312A		
Autoclave	CertoClav EL 18L GS 125/140°C	8501318		
Autoclave	CertoClav EL 18L GS 115/121°C	8501318A		
Accessories	Wire basket, stainless steel 23cm, with handle	8583010		
Accessories	Exhaust filter (0,01 Micron) complete with vessel	8582010		
Accessories	Condensation vessel (5L) with hose (1,0m)	8582020		
Accessories	Exhaust hose (1,0m)	8582016		
Accessories	Service Set 1	1250499		