



2mL crimp top vials




d1	11.6 ± 0.3
h1	32.0 ± 0.5
d2	10.7 ± 0.2
d4	6.0 ± 0.2
h4	3.6 ± 0.2
d3	8.3 ± 0.3
h3	≈ 7.5
S2	≥ 0.4
S1	0.85 ± 0.05

2mL screw top vials (8-425)



d1	11.6 ± 0.25
d2	≤ 8.0
d4	5.0 ± 0.30
h1	32.0 ± 0.50
d3	8.8 ± 0.30
h3	≈ 7.8
h5	1.6
h6	1.0
S2	≥ 0.50
S1	0.85 ± 0.04

2mL screw top vials (9-425)



d1	11.6 ± 0.2
d2	8.1-8.5
d4	6.0 ± 0.20
h1	32.0 ± 0.50
d3	9.4 ± 0.25
h3	7.3-7.9
h5	≈ 1.2
h6	≈ 0.9
S2	≥ 0.50
S1	0.85 ± 0.04