



Heating mantle | round flasks

for spherical flasks, integrated magnetic stirrer and controller, up to 450°C



1.000 ml round flask heating mantle

WHM for round flasks, 100°C-450°C with built-in magnetic stirrer integrated controller

Ideal for:

- spherical flasks 50 - 6,000 ml

Features:

- full choice of vessel, capacity, housing and control are offered
- nickel chrome heating element
- coil type: acid-resistant, long life time
- built-in temperature controller
- all kinds of round flasks can be used (Duran®/Kimax/Pyrex/Iwaki)
- CE certified and unique serial number for tracing

Controller:

- built-in **stirring and temperature controller**



2x stirrer bars included

Application (flask brand)	Flask capacity (ml)	Flask Ø (mm)	Watts (W)	Inside depth (mm)	Mantle outside Ø x H mm	Support rod Ø x H mm	Order number	
							with 230V	with 120V
common-use	50	49 - 52	80	40	Ø 145 x 140	Ø 13 x 440	DH.WHM12030	DH.WHM12040
common-use	100	63 - 66	80	40	Ø 145 x 140	Ø 13 x 440	DH.WHM12031	DH.WHM12041
common-use	250	83 - 86	140	50	Ø 145 x 140	Ø 13 x 440	DH.WHM12032	DH.WHM12042
common-use	500	102 - 105	240	65	Ø 190 x 165	Ø 13 x 440	DH.WHM12033	DH.WHM12043
common-use	1,000	130 - 133	350	75	Ø 190 x 190	Ø 13 x 440	DH.WHM12034	DH.WHM12044
for DU./Iw.	2,000	165 - 169	450	100	Ø 270 x 200	Ø 13 x 660	DH.WHM12035	DH.WHM12045
for Ki./Py.	2,000	161 - 164	450	100	Ø 270 x 200	Ø 13 x 660	DH.WHM120351	DH.WHM120451
for DU./Ki.	3,000	184 - 187	550	105	Ø 270 x 200	Ø 13 x 660	DH.WHM12036	DH.WHM12046
common-use	5,000	221 - 224	800	131	Ø 325 x 226	Ø 13 x 800	DH.WHM12037	DH.WHM12047
for DU./Ki.	6,000	234 - 237	800	131	Ø 325 x 226	Ø 13 x 800	DH.WHM12038	DH.WHM12048

Remarks: DU. = Duran®, Ki. = Kimax, Iw. = Iwaki, Py. = Pyrex

Accessories: Magnetic stirrer bars can be found on page 28