

# Chest freezer XLT C300



- Low cost of ownership
- Limited maintenance and on-site service
- Easily installable & movable

- Efficient cyclopentane insulation
- Eco-friendly, HCFC-free Nature R refrigerant
- ISO 9001, 14001 & 45001 Certified

## Specifications

### Technical

|                          |                  |
|--------------------------|------------------|
| Gross volume             | 296 litres       |
| Net Volume               | 284 litres       |
| Temp. range              | -40°C to -60°C   |
| Temp. probes             | 1/PT 1000        |
| Compressor               | Embraco          |
| Baskets (incl./optional) | 3/12             |
| Refrigerant              | Nature R         |
| Noise level              | 55 dB            |
| Power consumption 24/h   | 5,5 kWh          |
| Power consumption        | 2.007 kWh / year |
| Voltage                  | 230V             |
| Frequency                | 50 Hz            |

### Features

|                                          |          |
|------------------------------------------|----------|
| Audible and visual alarm                 | YES      |
| Castors                                  | 4        |
| Integrated lock                          | YES      |
| Inner lids                               | NO       |
| Defrost                                  | Manual   |
| Battery Backup                           | NO       |
| USB                                      | NO       |
| Data connection                          | NO       |
| Digital controller with display          | XR30CX   |
| Data logging                             | NO       |
| High/Low temperature alarm               | YES      |
| Porthole for external probe & CO2 backup | Ø12 mm   |
| Connection for remote alarm              | Optional |

### Material

|       |                             |
|-------|-----------------------------|
| Inner | Precoated white steel       |
| Outer | Galvanized pre-coated steel |

### Outer dimensions

|                                  |         |
|----------------------------------|---------|
| Height (with castors)            | 885 mm  |
| Width                            | 1260 mm |
| Depth (incl. condenser)          | 630 mm  |
| Depth (incl. condenser & handle) | 680 mm  |
| Net weight                       | 71 kg   |

### Inner dimensions

|                   |             |
|-------------------|-------------|
| Height            | 635/385 mm  |
| Width             | 1100/890 mm |
| Depth             | 440 mm      |
| Insulation        | 80 mm       |
| Insulation type   | PUR         |
| Internal humidity | <5 %        |

## XLT C300



Product no.

N112003

## Storage capacity

| Box height | No. of boxes | No. of racks | Total no. of boxes | Product no. |
|------------|--------------|--------------|--------------------|-------------|
| 50 mm      | 6            | 3            | 18                 | N311200     |
|            | 11           | 18           | 198                | N311220     |
|            |              |              | <b>= 216</b>       |             |
| 75 mm      | 4            | 3            | 12                 | N311300     |
|            | 8            | 18           | 144                | N311320     |
|            |              |              | <b>= 156</b>       |             |
| 100 mm     | 3            | 3            | 9                  | N311400     |
|            | 5            | 18           | 90                 | N311420     |
|            |              |              | <b>= 99</b>        |             |