



Universal Mini Centrifuge uniCFUGE 2



- Universal rotor for simultaneous use of centrifuge tubes and PCR-strips
- Two speeds selectable

Compact and reasonably priced mini centrifuge with **universal rotor** for 6 x 1.5/2.0 ml and 6 x 0.5 ml tubes and 2 times 8 x 0.2 ml PCR-strips.

Speed can be adjusted to 4000 rpm or 6000 rpm. Powerful motor for fast starting. Quick stop when the lid is opened.

Specifications

Max. RCF:	1200 xg and 2000 xg
Max. speed:	4000 rpm and 6000 rpm
Max. load:	6 x 0.5/1.5/2 ml and 2 x 8-tube PCR strips (0.2 ml)
Dimensions (diam. x H):	140 x 120 mm
Weight:	0.7 kg
Power supply:	100-240 V, 50/60 Hz, 1 A

Type	PK	Cat. No.
LLG-uniCFUGE 2 with universal rotor with EU-plug	1	6.263 510
LLG-uniCFUGE 2 with universal rotor with UK-plug	1	6.263 511



Compact Mini Centrifuge uniCFUGE 2/5



- Rotor for 4 x 5 ml tubes
- Two speeds selectable

Compact and reasonably priced mini centrifuge with **rotor for 4 x 5 ml tubes**.

Speed can be adjusted to 4000 rpm or 6000 rpm. Powerful motor for fast starting. Quick stop when the lid is opened.

Specifications

Max. RCF:	1200 xg and 2000 xg
Max. speed:	4000 rpm and 6000 rpm
Max. load:	4 x 5 ml tubes
Dimensions (diam. x H):	140 x 120 mm
Weight:	0.7 kg
Power supply:	100-240 V, 50/60 Hz, 1 A

Type	PK	Cat. No.
LLG-uniCFUGE 2/5 with rotor for 5ml Tubes with EU-plug	1	6.263 515
LLG-uniCFUGE 2/5 with rotor for 5ml Tubes with UK-plug	1	6.263 516