

Chemistry scoop



Scoop for aggressive liquids, complete with strong handle, entirely made of polypropylene.

- ▶ Rod length 100 cm, Ø 20 mm
- ▶ Without metal parts

Stainless steel scoop



In contrast to conventional scoops, the stainless steel scoop has no openings, undercuts or crevices in which residues could collect. Pore-free finishing produced by high-quality laser welding seams.

The stainless steel scoop can be flame-sterilised and can therefore be made completely germ-free.

- ▶ Scoop made of stainless steel V2A (1.4301), complete with rod
- ▶ Sterilisable by flame sterilisation or steam sterilisation

Chemistry scoop



Capacity ml	Price €	Quantity scale from	Scale price €	Item no.
250	55.60	10	50.20	5353-0250 N4
600	57.55	10	51.70	5353-0500 N4
1000	58.90	10	52.85	5353-1000 N4

Stainless steel scoop



Total length cm	Capacity ml	Price €	Item no.
115	1000	209.60	5354-1000 N4