




RED30 Series

RED30S





RED LINE - RED30 Series

RED30S

Optical System	Greenough
Observation Tube	Binocular head
Interpupillary distance	52-69mm
Diopter adjustment	On the left tube, +/- 5 diopter
Inclination	45° inclined
Eyepieces	WF10X/20mm
Objectives system	Turret changer with LED incident illumination
Objectives magnification	2X - 4X
Working distance	50mm
Stand	Fixed arm type stand
Base	225x179mm
Column/Arm	200mm high
Head holder	Fixed
Focus mechanism	Coarse focusing system with tension adjustment
Focusing stroke	47mm
Incident illumination	LED 3.5V/100mW; with intensity control
Transmitted illumination	LED 3.5V/100mW; with intensity control
Transformer	External
Power supply	Charger 100-240V or rechargeable batteries
Accessories included	Black, white and frosted stage plates, dust cover, charger and batteries
Length x Width x Height	225x179x250mm
Weight	3.5Kg