




**Special transfer siphon with two separate
LN2-fluid valves and two separate
ascending pipes for containers without
permanently installed safety fittings
(Type ALU / ALU-CD)**



EKI-2V



EKI-1V-1MAG



EKI-2MAG

Special transfer siphon with two separate LN2-fluid valves and two separate ascending pipes for containers without permanently installed safety fittings (Type ALU / ALU-CD)



EKI-2V
Art. No. 2611-D1

consisting of:

- 1 x EKI = transfer siphon with manometer and overpressure valve 0.5bar
- 2 x LN2-fluid valve with a screw adapter
- 1 x pressure-release valve

EKI-1V-1MAG
Art. No. 2611-D2

consisting of:

- 1 x EKI = transfer siphon with manometer and overpressure valve 0.5bar
- 1 x LN2-fluid valve with a screw adapter
- 1 x magnetic valve for LN2 with a screw adapter
- 1 x pressure-release valve



EKI-2MAG
Art. No. 2611-D3

consisting of:

- 1 x EKI = transfer siphon with manometer and overpressure valve 0.5bar
- 2 x magnetic valve for LN2 with a screw adapter
- 1 x pressure-release valve



Custom LN2 transfer siphons on request

We design and manufacture your custom-made LN2 transfer siphon according to your needs. Please do not hesitate to contact us: info@kgw-isotherm.de

