


Heating mantle | round flasks

for spherical flasks, external controller available, up to 450°C



1.000 ml round flask heating mantle

WHM for round flasks, 450°C for use with external controller

Ideal for:

- spherical flasks 50 - 6,000 ml
- use with magnetic stirrer

Features:

- full choice of vessel, capacity, housing and control are offered
- highly processed aluminium case offers robust standing, easy handling
- nickel chrome heating element
- coil type: acid-resistant, long life time
- K-type thermo sensor is basically embedded
- CE certified and unique serial number for tracing

Controller:

- optional remote thermo controllers are available (analog: WHM-C10A, digital: WHM-C10D)



WHM-C10D
digital controller



WHM-C10A
analog controller

Application (flask brand)	Flask capacity (ml)	Flask Ø (mm)	Watts (W)	Inside depth (mm)	Mantle outside Ø x H mm	Support rod Ø x H mm	Order number	
							with 230V	with 120V
common-use	50	49 - 52	80	48	Ø 145 x 75	Ø 13 x 440	DH.WHM12050	DH.WHM12060
common-use	100	63 - 66	80	50	Ø 145 x 75	Ø 13 x 440	DH.WHM12051	DH.WHM12061
common-use	250	83 - 86	140	50	Ø 145 x 75	Ø 13 x 440	DH.WHM12052	DH.WHM12062
common-use	500	102 - 105	240	85	Ø 190 x 110	Ø 13 x 440	DH.WHM12053	DH.WHM12063
common-use	1,000	130 - 133	350	85	Ø 190 x 110	Ø 13 x 440	DH.WHM12054	DH.WHM12064
for DU./Iw.	2,000	165 - 169	450	110	Ø 270 x 135	Ø 13 x 660	DH.WHM12055	DH.WHM12065
for Ki./Py.	2,000	161 - 164	450	110	Ø 270 x 135	Ø 13 x 660	DH.WHM120551	DH.WHM120651
for DU./Ki.	3,000	184 - 187	550	110	Ø 270 x 135	Ø 13 x 660	DH.WHM12056	DH.WHM12066
common-use	5,000	221 - 224	800	140	Ø 325 x 165	Ø 13 x 800	DH.WHM12057	DH.WHM12067
for DU./Ki.	6,000	234 - 237	800	140	Ø 325 x 165	Ø 13 x 800	DH.WHM12058	DH.WHM12068

Remarks: DU. = Duran®, Ki. = Kimax, Iw. = Iwaki, Py. = Pyrex

Accessories: Remote controller and temperature probes on page 46