MiniSampler | taking samples – comfortable, quick and clean sampling





With the vacuum sampler MiniSampler, the sample liquid only comes into contact with the suction hose and the original sample bottle. The possibility of cross-contamination is greatly reduced. Immediately after taking the sample, the bottles are sealed. The disposable suction hose can quickly be replaced. This guarantees that the samples are never contaminated. If necessary, a new hose can be used for every single sample.
Due to its small diameter (8 mm), the hose can reach narrow, seemingly inaccessible sampling areas.

Two models:


PE MiniSampler for general sampling, max. suction height 5 m*
PTFE MiniSampler for ultra-pure sampling, max. suction height 5 m*
* depending on the medium

- Also suitable for viscous oils
- Complete in a handy transport case

Use

1. Insert the hose and the bottle into the pump.
2. Insert hose end to desired sample depth.
3. Move the pump piston back and forth and collect the sample. The liquid flows directly into the bottle.
4. Take off the filled bottle, close and label it.

⚠ Not for aggressive acids. Vapors could destroy pump.



MiniSampler PE | for general sampling



Ideal for general sampling from barrels, canisters, tanks etc.
MiniSampler PE complete with accessories in transport case, contents: Vacuum pump, 10 PE bottles, 100 ml, 10 m PE suction hose 6 x 8 mm, 1 hose cutter, 1 stainless steel hose weight made of V2A (1.4301), 20 seal-it bottle seals.

MiniSampler PE

Price €	Quantity scale from	Scale price €	Item no.	
154.30	5	137.00	5305-0100	N4

Accessories

Description	Price €	Item no.	
MiniSampler pump with adapter	66.85	5305-0010	N4
MiniSampler standard pump	42.60	5305-0011	N4
MiniSampler adapter for glass thread GL 45	118.60	5305-0209	N4
MiniSampler adapter for sample tubes PFA 6205-0090, 6202-0180	72.50	5305-1020	N2
Hose cutter	7.00	0006-9506	N3
Sampling bottle HDPE, 100 ml	1.50	0340-0100	N3
Sampling bottle HDPE, 250 ml	1.83	0340-0250	N3
Sampling bottle HDPE, 500 ml	2.15	0340-0500	N3
Sampling bottle HDPE, 1000 ml	2.78	0340-1000	N3
Suction hose PE on roll 6 x 8 mm, roll with 10 m	13.35	8805-0608	N3
Suction hose PE on roll 6 x 8 mm, roll with 100 m	41.95	8878-0608	N3

MiniSampler PTFE | for ultrapure sampling



Ideal when maximum purity is needed.
MiniSampler PTFE complete with accessories in its transport case, contents: Vacuum pump, 1 sample bottle PFA 180 ml, 5 m PTFE suction hose 6 x 8 mm, 1 hose cutter, 1 stainless steel hose weight made of V2A (1.4301), 20 seal-it bottle seals.

MiniSampler PTFE

Price €	Quantity scale from	Scale price €	Item no.	
253.35	5	230.20	5305-0110	N4

Accessories

Description	Price €	Item no.	
MiniSampler pump with adapter	66.85	5305-0010	N4
MiniSampler standard pump	42.60	5305-0011	N4
MiniSampler adapter for glass thread GL 45	118.60	5305-0209	N4
MiniSampler adapter for sample tubes PFA 6205-0090, 6202-0180	72.50	5305-1020	N2
Hose cutter	7.00	0006-9506	N3
Sample boxes PFA, 90 ml	52.60	6205-0090	N3
Sample boxes PFA, 180 ml	68.10	6205-0180	N3
Suction hose PTFE on roll 6 x 8 mm, roll with 5 m	37.65	8804-0608	N3