1



1 LLG-Exact-Temp thermometers

For monitoring critical temperature conditions in
- deep freezes
- refrigerators
- cold rooms
- blood bank refrigerators
- incubators
- warming chambers

Characteristics:
- thermometer certifiable to National Standards
- supplied with works test certificate
- thermometer stem safely encapsulated in PTFE
- each thermometer is mounted in a square-shaped shatter-proof plastic bottle filled with a liquid to ensure accurate readings despite the opening and closing of instrument doors.
- each thermometer has an individual serial number.

LLG-Exact-Temp Standard thermometers, red spirit-filled

For	Measuring range °C	Grad. °C	Length mm	Bottle filling	PK	Cat. No.
Freezers	-90 ... +20	1.0	145	Sand	1	**9.235 663**
Refrigerators	-30 ... 0	0.5	125	Ethylene Glycol	1	**9.235 664**
Blood banks	-5 ... +20	0.5	145	Ethylene Glycol	1	**9.235 665**
Coolers	-5 ... +15	0.5	125	Ethylene Glycol	1	**9.235 666**
Incubators	+18 ... +50	0.5	135	Ethylene Glycol	1	**9.235 667**
Warming rooms, ovens	+20 ... +130	1.0	135	Vermiculite	1	**9.235 668**
Warming rooms, ovens	+35 ... +200	1.0	145	Vermiculite	1	**9.235 669**

LLG-Exact-Temp Standard thermometers, mercury-filled

For	Measuring range °C	Grad. °C	Length mm	Bottle filling	PK	Cat. No.
Freezers	-30 ... 0	0.5	115	Ethylene Glycol	1	**9.235 670**
Refrigerators	-5 ... +15	0.5	125	Ethylene Glycol	1	**9.235 671**
Incubators	+18 ... +50	0.5	135	Ethylene Glycol	1	**9.235 672**
Warming rooms, ovens	+20 ... +130	1.0	135	Vermiculite	1	**9.235 673**
Warming rooms, ovens	+35 ... +200	1.0	145	Vermiculite	1	**9.235 674**

LLG-Precision thermometers Exact-Temp, blue spirit-filled

For	Measuring range °C	Grad. °C	Length mm	Bottle filling	PK	Cat. No.
Freezers	-25 ... -5	0.1	210	Ethylene Glycol	1	**9.235 675**
Refrigerators	-2 ... +10	0.1	180	Ethylene Glycol	1	**9.235 676**
Blood banks	-2 ... +10	0.1	180	Ethylene Glycol	1	**9.235 677**
Incubators	+15 ... +30	0.1	180	Ethylene Glycol	1	**9.235 678**
Incubators	+25 ... +45	0.1	190	Ethylene Glycol	1	**9.235 679**
Warming rooms, ovens	+50 ... +115	0.1	200	Vermiculite	1	**9.235 680**

LLG-Precision thermometers Exact-Temp, mercury-filled

For	Measuring range °C	Grad. °C	Length mm	Bottle filling	PK	Cat. No.
Cold rooms, Freezers	-25 ... -5	0.1	210	Ethylene Glycol	1	**9.235 681**
Refrigerators	-2 ... +10	0.1	180	Ethylene Glycol	1	**9.235 682**
Blood banks	-2 ... +10	0.1	180	Ethylene Glycol	1	**9.235 683**
Incubators	+15 ... +30	0.1	180	Ethylene Glycol	1	**9.235 684**
Incubators	+25 ... +45	0.1	190	Ethylene Glycol	1	**9.235 685**
Warming rooms, ovens	+95 ... +115	0.1	200	Vermiculite	1	**9.235 686**